ADVANCED WOUND CARE IN AFRICA: TREATMENT OF CHRONIC, INFECTED FOOT ULCER WITH MELMAX.

t=0



May 2016
Bedridden female
patient (age ± 55) has
been struggling with a
stagnating wound since
1994.

Several treatment modalities have been tried, all unsuccessful. Amputation is expected.

Provided by:



 $t = \pm 45 \text{ days}$



July 2016 Application of MelMax; dressing change three times weekly; rinse with NaCl twice weekly.



 $t = \pm 90 \text{ days}$



August 2016
Superabsorbent
secondary dressing
used, fixated with
regular elastic gauze.

Satisfactory drainage result; tissue returning to normal (from redness, pain, swollen).

Granulation tissue and first re-epithelialisation tissue visible.

 $t = \pm 120 \text{ days}$



September 2016
Complete
epithelialisation. A
hydrocolloid patch was
applied temporarily to
prevent further
mechanical trauma.
Patient was able to walk
(with support) for the
first time in years.

