Principelle IF

HONEY BASED IONOGEN FORMULATION OINTMENT

Doc. No.1608-01.02

UK –Principelle IF® Tube

Product Description

Principelle IF is a sterile ointment containing a blend of trace elements and honey. Minerals, trace elements and oxides are directly blended with medical grade dark buckwheat honey. The formulation is gradually absorbed by the wound surface after application. Principelle IF® ointment can be used for all types of wounds, including burns. Tube packaging especially suitable for use in cavity wounds.

Product Characteristics

Supervised by a healthcare professional, Principelle IF® is indicated for the treatment of acute and chronic wounds, including contaminated and infected wounds. It may be used on moderately and lightly exuding wounds, and dry wounds. The wound needs to be covered with a secondary dressing. Depending on the level of exudate production an absorptive dressing is recommended. This can be a traditional absorptive dressing or Principelle Matrix®. Strongly exuding wounds reduce wear time. Principelle IF® normalises the wound micro-environment by restoring the protease balance and regulation of bacterial contamination in therapy resistant wounds. This results in better conditions for re-epithelialisation. Proteases play an important role in tissue degradation and regeneration and in the development of epithelium. Modulation of the wound micro-environment is associated with appropriate stages of healing. The characteristics of Principelle IF® are also directed at the regulation of moisture content in the micro-environment of the wound and at the regulation of the pH on the wound surface during the healing process. Principelle IF® does not inhibit the drainage of wound exudates.

Indications for Principelle IF® tube

Acute wounds:

- Burns
- Surgical wounds
- Traumatic wounds

Other acute wounds where fast epithelialisation is required

Chronic wounds:

- Decubitus
- Cavities
- Ulcera cruris
- Diabetic ulcers

In conjunction with corticosteroid use

Contra-indications

Principelle IF® may not be used on patients with a known hypersensitivity to Ionogen Formulation or honey.

Instructions for Use

Application:

1. Carefully inspect the wound.
2. Cleanse the wound according to standard instructions.
3. Open the tube.

If the tube has been opened or damaged, then consider the Principelle IF® tube to be non-sterile.

4. Gently apply Principelle IF® and ensure that this leads to maximum contact with the wound.

Note:

In cases of moist or wet wounds, Principelle IF® can be covered with a secondary dressing, which will result in the retention of wound exudate.

For Single Use Only

Principelle IF® is sterilised with gamma radiation. The sterility of Principelle IF® is only guaranteed for unopened and undamaged packages. Principelle IF® is not to be re-sterilised.

Principelle IF® tube Packaging	Contents	Order number
1 tube	10 g	PPIF-10
1 tube	20 g	PPIF-20

Storage:

Store at room temperature (5-25 °C), dry, and keep away from bright light.

	For Single Use
	Sterile (Gamma)
	Store in a dry place
	Do not store in direct sunlight

DO NOT use if tube is damaged or opened prior to use.

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