

Your answer to colonised and infected wounds

MelMax® cleaning and healing wounds

MelMax® is a sterile, inert wound dressing impregnated with a mixture of MMP-regulating metal-ion formulation and buckwheat honey specially selected for medical use.


Metal-ion formulation, a blend of trace elements that normally occur in wound exudates and serum, aids regulating MMP production in the wound bed.

Buckwheat honey is honey with a rich source of phenolic compounds shown to be responsible for superior anti-oxidant activity ¹⁾.

Combining buckwheat honey with the metalion formulation is the answer to colonised and infected wounds with diminished healing capacity.

Since honey will clean colonised and infected wounds and the metal-ion formulation facilitates wound closure, MelMax® is the product of choice.





MelMax Me

